

Water determination in liquids and gases

Universal basic instrument
for coulometric Karl Fischer titration



aqua 40.00 **Basic Module**

made by **ECHⁱ**

aqua 40.00 Basic Module

Water determination in liquids and gases

Product description

The **AQUA 40.00 Basic Module** determines the water content of liquid samples and gases quickly and precisely. The instrument is based on coulometric Karl Fischer titration and requires little space in the laboratory due to its compact design.

The sample is simply dispensed through a septum directly into the measuring cell. The cell is optimally designed and exhibits extremely low background drift. The titrator is thus particularly suitable for water determination in the trace range. The titration speed is automatically adjusted to the actual amount of water to be titrated because of continuously controlled electrolysis current. This results in very short analysis times.

All common Karl Fischer reagents can be used for the titrator. In many applications,

a conventional diaphragm is not required in the generator electrode. Thus, only one coulometric reagent is required.

The application range of the **AQUA 40.00 Basic Module** can be extended with various additional modules. These modules can be retrofitted. This makes the analyzer system suitable for liquid and gaseous samples as well as for solid and pasty samples.



AQUA 40.00 Basic Module

Advantages

- No sample preparation
- Dosing by syringe directly into the titration cell
- Very short analysis times
- Suitable for all common Karl Fischer reagents
- Basic Module with low space requirement
- Additional modules available for different sample types

Applications

The AQUA 40.00 Basic Module is suitable for samples that can be dosed directly into the titration cell, e. g.

- Solvents
- Oils and fats
- Fuels, biodiesel
- Hydraulic and brake fluids
- Gases

Features

- Fast and reproducible analyses
- Continuously controlled electrolysis current
- Extremely low background drift
- Simple, clear software
- User-friendly handling
- Status display through large colouring


Evaluation of titration test according to ISO 9001-9003

No.	Sample amount [mg]	Dosed quantity [µg]	Detect. quantity [µg]	Detect. Water
1	459.85	459	461	1.003
2	425.45	425	422	0.991
3	505.20	504	505	1.000

Given set value:	0.998	µg/mg
Detect. actual value (mean):	0.999	µg/mg
Detect. reproducibility:	Rel. standard tolerance:	0.53 %
Detect. accuracy:	Rel. accuracy:	0.11 %
Test criterion met:	Yes	

Close Help

Evaluation of the titration test



Typical measurement of a water standard

Reagent consumption

General | Electrolyte

Parameter

Name: Hydranal-Coulomat AG

Charge: 5243C

Volume: 100.0 ml

Capacity: 10.0 mg/ml

Reagent consumption

Completed meas.: 7

Duration: 105 days

Conversion: 30.417 mg

OK Cancel Help

Reagent consumption

Additional modules

The AQUA 40.00 Basic Module can be expanded to the AQUA 40.00 Vario. This allows the water determination in all sample types by using headspace technology.

The AQUA 40.00 Vario is available as a manual version or with an autosampler. Both versions are suitable for headspace vial sizes 2 R - 50 R and can be easily adapted to the respective vial size.

If high heating temperatures (up to 1300 °C) are required, the AQUA 40.00 Basic Module can be combined with the high temperature oven HT 1300. Then samples such as inorganic salts, building materials, metals, molecular sieves, oxides/hydroxides can also be analysed.



AQUA 40.00 Basic Module with connection module and high temperature oven HT 1300



AQUA 40.00 Basic Module extended to AQUA 40.00 Vario as manual version



AQUA 40.00 Basic Module extended to AQUA 40.00 Vario PLUS, the automatic version with sampler

Accessories

The **SWOP BOX** makes it easy and convenient to change reagents in titration cells. The module can be used, for example, on any Karl Fischer titrator of any design if the titration cell has at least one free connection.



Universal reagent exchange module SWOP Box

An automatic Liquid sampler is available for the AQUA 40.00 Basic Module. It is controlled by the software of the AQUA 40.00 Basic Module. Large sample quantities can be measured much faster with it.



Liquid sampler as accessory for the AQUA 40.00 Basic Module

Technical specifications of the AQUA 40.00 Basic Module

Measuring range: 1 µg ... 100 mg absolute
 Resolution: 0.1 µg
 Reproducibility: ± 3 µg for 10 ... 1000 µg, 3 % for > 1 mg
 Generator current: 0 ... 250 mA
 Sample volume: 0.01 ... 20 mL (direct injection)
 Volume of reagent: 100 mL
 Typical measuring time: 5 min (depending on the water content)
 Result display: µg, µg/L, mg/L, ppm, %, mC, customized with formula generator

Indication: biamprometrical, polarization with square-wave voltage
 Power supply: 230 V/50 Hz or 115 V/60 Hz
 Interface: RS 232
 Balance connection: RS 232
 Dimensions: 112 x 448 x 208 mm (W x H x D)
 Weight: 5 kg
 Device control: PC software (PC not included in the scope of delivery)

Book your online demo in the ECH Studio

ECH Scientific have a state of the art laboratory fitted with online presentation capabilities, allowing us to bring product demonstrations live and in full HD, with multiple camera angles and software sharing capabilities enabling us to deliver a full demo experience remotely. Please contact info@echscientific.com to book your session.

ECH Elektrochemie Halle GmbH

Otto-Eißfeldt-Str. 8
 D-06120 Halle (Saale)
 Germany

Tel.: +49 (0) 345 279570-0
 Fax: +49 (0) 345 279570-99

Email: info@echscientific.com • www.ech.de • www.aquamaxkf.com

ECH Scientific Limited

Building 69, Wrest Park, Silsoe
 Bedfordshire, MK45 4HS
 United Kingdom

Tel.: +44 (0) 1525 404747
 Fax: +44 (0) 1525 404848

ECHⁱ
SCIENTIFIC
 part of ECH Elektrochemie Halle
 Global Sales Division

the ECH advantage

in-lab | mobile | on-line | process