

Manual Oil Sampling Kit

made by **ECH**

Product description

During the sampling of oils you must avoid any contact with ambient air.

The Manual Oil Sampling Kit is designed for contamination-free taking of oil samples.

A constant temperature of oil sample is important for correct analysis. The temperature is automatically measured during flow of oil through the Manual Oil Sampling Kit. When the temperature of oil sample is constant you can fill the sample into the gas tight glass syringe.

The Manual Oil Sampling Kit is easy to use. Versions for various syringe sizes are available (10 mL, 50 mL, 100 mL).

Applications

- Mineral oils
- Brake fluids
- Biological oils
- Food- and pharmaceutical products
- Vegetable oils
- Oxidation sensitive oils

Features

- Gas tight sampling kit for oils
- No contact of oil sample with ambient air
- Versions for various syringe sizes
(10 mL, 50 mL, 100 mL)



Manual Oil Sampling Kit with glass syringe (100 mL)